



Plant Biotechnology

Bayer Russia Plant Biotechnology
Conference

July 2023 – week 1





Bayer Russia Plant Biotechnology Conference:

Day 1	Plant Biotechnology Overview and Evolution
Day 2	Introduction to Plant Vectors & Agrobacterium
Day 3	Introduction to Plant Transformation - Corn
Day 4	Soy Plant Transformation & Gene Expression
Day 5	Introduction to Plant Health & Controlled Environment (CE)



Introduction to vectors

Xudong Ye, PhD.
Bayer Crop Science



About me:

Principal Scientist, Bayer Science Fellow

PhD Institute of Plant Science, ETHZ, Switzerland (Prof. Ingo Potrykus)

Ryegrass transformation, fructan metabolism

M.S. Plant Pathology, Nanjing Agricultural University, China

B.S Plant Protection, Zhejiang (Agricultural) University, China

My research fields

Plant transformation (dicot and monocot, nuclear and plastid)

Agrobacterium biology (strain modification, vector strategies)

Molecular biology

Gene editing

Patents


Publications



Agenda

- Vectors designs for particle gun and *Agrobacterium*
- Selectable markers
- Scorable marker genes (GFP, GUS)
- Backbone
- Main cloning strategies

Plant gene structure




From genome to plasmid: cloning DNA fragments by in vitro manipulation




Plasmid

Vector

Construct



Plasmids for gun bombardment are different from *Agrobacterium* vectors



- Smaller than *Agrobacterium* binary vectors
- Single origin of replication, high copy for higher yield
- May need to release the GOI by restriction enzymes and HPLC purify the gene fragment for product development





Binary vectors for *Agrobacterium*-mediated transformation

Free vector in Agrobacterium, replicate in E. coli and Agrobacterium

Necessary for T-DNA delivery

VirE1-3

VirD1-5


VirC1-2

VirG (activator)


VirB1-11

VirA (sense phenolic or acidic environment)


pTiC58




A typical binary vector for *Agrobacterium*-mediated transformation



Destructive assay



GFP transient expression



[Plant development inhibitory genes in binary vector backbone improve quality event efficiency in soybean transformation | SpringerLink](#)

[Enhanced production of single copy backbone-free transgenic plants in multiple crop species using binary vectors with a pRi replication origin in *Agrobacterium tumefaciens* | SpringerLink](#)

Frequently used selectable marker genes

Selectable marker gene	Selection reagent	Comments
<i>nptII</i>	Kanamycin, G418, paromomycin	Antibiotic, lethal
<i>epsps-CP4</i>	glyphosate	Herbicide, lethal
<i>bar, pat</i>	Glufosinate, phosphinothricin	Herbicide, lethal
<i>hpt</i>	hygromycin	Antibiotic, lethal
<i>PMI (manA)</i>	mannose	Non-lethal
<i>aadA</i>	Spectinomycin, streptomycin	Antibiotic, non-lethal



Introduction to cloning methods

The bible of cloning:

Molecular Cloning: A Laboratory Manual (Fourth Edition) (cshlpress.com)

- 1) Traditional cloning method by ligation of compatible fragments with T4 DNA ligase
- 2) Golden-gate cloning for multiple fragment ligation using type II restriction enzymes
- 3) Assemble multiple fragments by fusion with short homologue ends (Gibson assembly, hot fusion, In-fusion etc.)



Restriction enzymes

<https://www.pnas.org/content/102/17/5905>

How restriction enzymes became the workhorses of molecular biology

<https://www.nature.com/scitable/topicpage/restriction-enzymes-545/>

Restriction Enzymes

<https://www.britannica.com/science/restriction-enzyme>

[Restriction enzymes \(ppt\)](#)

<https://capricorn.bc.edu › 11-Restriction-mapping>

5'	GAATTC	3'	
3'	CTTAAG	5'	<i>EcoRI</i>

5'	GGTACC	3'	
3'	CCATGG	5'	<i>KpnI</i>

5'	AGGCCT	3'	
3'	TCCGGA	5'	<i>StuI</i>

Purchase restriction enzymes

www.NEB.com

<https://www.thermofisher.com/>

Different REs may have different requirement:

- Buffer (salts, pH)
- Temperature (37, 25, 55, 65 °C)
- Double digestion consideration
- DNA methylation

Dam (GATC) (e.g. XbaI)

Dcm (CCAGG and CCTGG)

Dam- E coli is available for Dam- DNA (SCS110, INV110)

<https://blog.addgene.org/plasmids-101-methylation-and-restriction-enzymes>

Read user manual for instruction


No double dip when take RE

Cloning method 1: restriction enzyme digestion and ligation

Considerations for ligation-based cloning



- Compatible cut ends between fragments, if can't find, need to **blunt** the fragment ends by T4 DNA polymerase
- Prevent vector backbone self-ligation
 - if vector cut ends are compatible, **dephosphorylate** the vector fragment (e.g. CIP, or Shrimp phosphatase)
 - select 2 incompatible restriction cut for cloning to avoid
- Always purify fragments by electrophoresis
- Adjust vector fragment/insert ration to 1:2 or 1:3**
- DNA loads: for 15-25 kb vector backbone, use 50~100 ng/20 µl reaction volume, adjust insert in proportion
 - 3~15 kb 5~50 µg/20 µl reaction volume
 - Blunt end ligations use higher end concentration
 - Stick end ligations use lower end concentration
- Ligation at 16 °C overnight is recommended




<i>EcoRI</i>	<i>EcoRI</i>	<i>EcoRI</i>	<i>EcoRI</i>	<i>KpnI</i>	<i>KpnI</i>
5' GAATTC 3'	5' p-AATTC — G-OH 3'	5' p-AATTC — GGTAC-OH 3'	5' p-AATTC — GGTAC-OH 3'	5' GGTACC 3'	
3' CTTAAG 5'	3' HO-G — CTTA-p 5'	3' HO-G — C-p 5'	3' HO-G — C-p 5'	3' CCATGG 5'	

T4 DNA ligase ligates compatible DNA end together



<https://www.neb.com/products/m0202-t4-dna-ligase#Product%20Information>

How to calculate the molarity of DNA ends



- 1) Digest 2 µg pUC18 with SphI/EcoRI in 50 µl volume
- 2) Digest the insert DNA template with the same enzymes (2~5 µg/50 µl)
- 3) Fractionate digested plasmids in 1% Agarose gel
- 4) Excise the expected size bands with a scalpel over UV light box or Darkbox
- 5) Use commercial kit to dissolve the Agarose and purify the DNA fragments
- 6) Elute the DNA with water
- 7) Measure the DNA concentration
- 8) Prepare ligation reaction:
 - 1) 2 µl 10x T4 ligase buffer
 - 2) 5 ng pUC18 *EcoRI* / *SphI*
 - 3) 2x or 3x end moles of the insert
 - 4) Add H₂O to final 19 µl
 - 5) Add 1 µl T4 DNA ligase
 - 6) Mix well with pipet, incubate at 16 °C overnight





Transfer the ligate into E coli competent cells

https://blog.addgene.org/plasmids-101-common-lab-e-coli-strains?gclid=EAlaIQobChMlu_-bkeDi7gIVCjiGCh0ILA-nEAAYASAAEgLy2_D_BwE

DH10B is the most common strain used by us

ElectroMAX™ DH10B Cells

<https://www.thermofisher.com/order/catalog/product/18290015#/18290015>



Electroporation of plasmid into the competent cells

1 ul ligate into 20 ul aliquot of ice-thawed competent cell in an Eppendorf tube

Transfer into 1 mm gap BioRad cuvette

Tap on surface to be sure the drop is spread in the cuvette (prevent spark)

Set up 1.8 kv




Press pulse button

Add 500-1000 ul SOC, transfer into 14-ml tube

shake for 45-1 hr

Take 10 ul and 200 ul and spread into 2 plates with corresponding plasmid antibiotic


Reverse plates, culture at 37 °C overnight





Plasmid miniprep (alkaline lysis)

- Aliquot **2** ml LB with corresponding antibiotic (kan 50 or spec 50 mg/L) for same day prep at 8~9 a.m., or **3** ~5 ml at 4~5 p.m. for overnight culture, in 14 ml round bottom tubes
- Pick a well separated colony with a toothpick tip, throw it into media, **8~12** individual colonies/construct
 - Shake at 37 °C for 5-6 hrs (2 ml same day) or overnight (3~5 ml)
 - Pour ~1.5 ml into an Eppendorf tube carefully, closed
 - Spin 10~30 seconds at maximum speed
 - Remove supernatant with a tip connected to a vacuum flask
 - Add **200 ul P1** (with RNase added, Qiagen maxipep kit) or individual order. Vortex to suspend
 - Add **200 ul P2** to lysate bacterial cells
 - Add **200 ul P3** and close lid, sandwich with racks and shake to neutralize the lysates
 - Spin **1**~3 min at maximal speed
 - Transfer the supernatants into fresh Eppendorf tubes (pre-filled with 600 ul isopropanol)
 - Close lids, sandwich with racks, shake to mix well
 - Spin for **5**~10 min at maximum speed
 - Remove the supernatants carefully with vacuumed tip without touching the pellets
 - Add 1.5 ml 70% ethanol, remove with tip right away to reduce salt residues
 - Dried by SpeedVac for 1-3 min, or leave open at room temperature for 30 min or overnight
 - Resuspend in **50~100** ul pure water.



The quality is good for subcloning or direct sequencing


Miniprep kits available from Qiagen or Promega
QIAprep Spin Miniprep Kit



Cloning method 2: Golden Gate cloning



Type IIS restriction enzymes are unique from "traditional" restriction enzymes in that they cleave outside of their recognition sequence, creating four base flanking overhangs

BbsI gaagac123456
(BpiI) cttctg123456



BsaI
5' - GGTCTCN
3' - CCAGAGN NNNN - 5'

SapI
5' - GCTCTTCN
3' - CGAGAAGN NNN - 5'





GGTCTCN **GGAG** AATGNGAGACC
CCAGAGNCCTC **TTAC** NCTCTGG

GGTCTCN **AATG** GCTTNGAGACC
CCAGAGNTTAC **CGAA** NCTCTGG

GGAG **AATG** **GCTT**
TTAC **CGAA** **TCAG**

GGAG **AATG** **GCTT**
TTAC **CGAA** **TCAG**

Cloning method 3: Gibson assembly or hot fusion





Gibson assembly <https://doi.org/10.1126/science.1190719>

Hot fusion <https://doi.org/10.1371/journal.pone.0115318>



Primer design example for Gibson assembly or hot fusion



Use 1~100 ng template DNA, amplify the 3.8 kb fragment for 25~30 cycles, gel purified. Mix 2 fragment at 1:2~3 ratio for fusion



Cloning by hot fusion method

<https://journals.plos.org/plosone/article?id=10.1371/journal.pone.0115318>

- Open recipient vector with restriction enzymes
- Prepare insertion fragment with PCR amplification which should have 20-25 bp overlap sequence in the vector
- Measure DNA concentration by Nanodrop
- Calculate DNA end mole (calculator is available from NEB or Promega)
- **Vector DNA 30-60 ng (~20 kb) + insert = 1:2 OR 1:3 end concentration in 10 ul**
- Add 10 ul **2x** hot fusion buffer, mix
- Place on PCR machine to complete the reaction cycles (50 °C for 1 hr, 20 °C for 10 min, 4 °C forever)
- Transfer 1 µl into 20 µl *E coli* competent cells on ice
- Electroporation
- Shake 1 hr
- Plate
- 37 °C overnight
- Pick single colonies
- Miniprep
- Sequencing verify

Example to calculate end molarity

DNA	Size (kb)	Weight (ng)
Vector	15	60
Fragment 1	6	72
Fragment 2	1	12

required mass insert (g) = desired insert/vector molar ratio x mass of vector (g) x ratio of insert to vector lengths

How to calculate the molarity of ends

<https://nebiocalculator.neb.com/#!/ligation>

<https://www.promega.com/resources/tools/biomath/>

$6/15 \times 60 \times 3 = 72$ ng (1:3 ratio)

$1/15 \times 60 \times 3 = 12$ ng (1:3 ratio)



A typical flowchart for vector construction

- Design vector map
- Request base vector and templates
- Primer design and order
- Digest base vector, gel purify, measure concentration
- PCR insert fragments or restrict cut out fragments to be cloned
- Gel purify fragments, measure concentration
- Hot fusion PCR fragments with the base vector fragment
- Electroporate 1 μ l to *E. coli* competent cells
- Plate overnight
- Pick 8~12 single colonies, miniprep
- Sequence verification
- Re-prep a verified clone/construct
- Archive the plasmid
- Register workplan for plant transformation
- Design molecular assay (copy #, backbone, edits, etc.)
- Data collection




Protocol training:

*Cpf1 fragment cloning into
an Agrobacterium binary vector*

Xudong Ye

Corn and soybean base vector for *LbCpf1* cloning



Cpf1 is synthesized



How to calculate the molarity of DNA ends

The best ratio for insert 1 fragment in a vector may be 1:2 or 1:3

Ligation Calculator

This tool will calculate the mass of insert required at several molar insert:vector ratios in the range needed for typical ligation reactions.

Ligation

Tutorials

Insert DNA length

4

kb

Vector DNA length

15

kb

Vector DNA mass

50

ng

Required insert DNA mass

13.33 ng (1:1)



26.67 ng (2:1)



40.00 ng (3:1)

66.67 ng (5:1)

93.33 ng (7:1)

Vector SpeI/StuI: ~15 kb

Cpf1 fragment: ~4 kb

Recommend using 50~100 ng/20 μ l reaction volume of vector backbone fragment when cloning a 15~25 kb vector



Digestion and ligation mixtures

SpeI/StuI	Backbone size	Cpf1 size	
Base vector	14.3 kb	--	
Cpf1 plasmid	14.3 kb	3.75kb	


SpeI/StuI 1 hr 37 °C, run gel

4 µg/50 µl **Base vector**
5 µg/50 µl **Cpf1** plasmid

Yield: **30 ng/µl** in 75 µl = 2.25 µg (recovery rate = 56% (2.25 µl / 4 µg))

Backbone fragment 25 ng/ µl in 75 µl = 1.87 µg

Cpf1 fragment **6.8 ng/ µl** in 75 µl = **0.51 µg**



[How to calculate the molarity of ends](#)

(Cpf1 fragment **size** proportion: 3.75 kb/14.3 kb=0.26)

Example: 1:2 ratio: 60 ng vector backbone : **31.2 ng** (60x0.26 (**3.75/14.3 kb**) x2) Cpf1 fragment

1:3 ratio: 60 ng vector backbone : **46.8 ng** (60x0.26 x3) Cpf1 fragment

Ligation in a 1.5 ml Eppendorf tube or 200 µl PCR tube:


- 1) 2 µl (x **30 ng**) **Base vector** SpeI/StuI in **20 µl**, control
- 2) 5 µl (x **6.8 ng**) **Cpf1 fragment** SpeI/StuI (~**34 ng**) in 20 µl (~1:2) **OR**
- 3) 7 µl (x **6.8 ng**) **Cpf1 fragment** SpeI/StuI (~**48 ng**) in 20 µl (~1:3)
- 4) 2 µl 10x T4 DNA ligase buffer
- 5) 1 µl 10x T4 DNA ligase **water to final 20 µl**

Ligate overnight, in a centrifuge setting at 16°C **OR** use PCR tube for reaction overnight in a PCR machine

1 µl is transferred into 20 µl DH10B cells with electroporation, add 1 ml SOC, shake **1** hr at 37°C

plate **20** and **200** µl into 2 plates, respectively.

Cut gel, purify with Promega column





Forgot in shaker for **>3 hours** before plating, therefore, plate less volume to avoid crowdy colonies

Upper plates: 100 ul
Lower plates: 10 ul





Longer digestion or CIP treat backbone fragment may decrease background but not necessary



60 ng backbone

60ng+34 ng Cpf1 (1:2)

60 ng+48 ng Cpf1 (1:3)




~72


Checking 8 colonies for Cpf1 inserts with *Stu*I/*Spe*I digestion




8 colonies in 2 ml LB with spec 75 mg/L, shook for 6 hrs, pour into 1.5 ml tube, miniprep with 200 μ l P1+P2+P3, wash with 70% EtOH, dried and suspended in 100 μ l water. Take 10 μ l suspension, digest with *Stu*I/*Spe*I in 20 μ l final volume for 30 min



Mixture method



Invert the sandwiched plates 3~6 times before centrifuge



Good for suspending pellets

Often too much for mixing, shearing DNA



What's next

- Bulk up plasmid DNA by inoculating the verified clones in large volume of LB medium with 75 mg/L spec selection
- Purify Plasmid DNA using commercial kit
- Suspend in water or 1x TE buffer, measure DNA concentration
- **Submit samples for sequence verification**
- Transfer the plasmid into *Agrobacterium* competent cells for plant transformation



Introduction to Agrobacterium



Outline

Agrobacterium tumefaciens and crown gall disease

- Discover bacterium pathogenicity

- Bacteria-free tumors

- Ti plasmid discovery

- Virulence gene induction

- The mechanism: Type IV secretion system (T4SS) and T-strand transfer

Agrobacterium as a vector for plant transformation

- Breakthroughs in plant genetic engineering

- Disarmed *Agrobacterium* strains


- Co-integration vector, binary vector, 2T-DNA vector

Discover crown gall disease and Agrobacterium

- The first written record of crown gall disease, on grape, dates from **1853**
Over-growths at wound sites and severely limits crop yields and growth vigor
- Fridiano Cavara (**1897**) found that a bacterium causes crown gall in grape
- **1907** - Erwin Smith and C.O. Townsend isolated a bacterium from galls on daisy. When inoculated the bacterium onto healthy plants, galls were reproduced





Erwin F. Smith






Crown gall production by *bacteria-free* tumor tissues

White, P.R. and Braun, A.C. (1941). Crown gall production by bacteria-free tumor tissues. Science. 94: [239-241](#);




Braun made fundamental discoveries about how *Agrobacterium* transforms plant cells

A few days after inoculation, tumors become independent of *Agrobacterium*




Conclusion: Viable bacteria are no longer necessary beyond two days post-inoculation. After this period, tumors become independent of the bacteria, because the bacteria have altered the host cells, by transferring some **factors** into them.

Virulence can be transferred from one *Agrobacterium* to another



A large plasmid in gall-inducing *Agrobacterium* transfers virulence in 1974



Zaenen, I., van Larebeke, N., Teuchy, H., **van Montagu, M. and Schell, J.** (1974). **Supercoiled circular DNA in crown-gall inducing *Agrobacterium* strains.** Journal of Molecular Biology. 86: [109-127](#).

Larebeke, N.V., Engler, G., Holsters, M., Den Elsacker, S.V., Zaenen, I., **Schilperoort, R.A.** and Schell, J. (1974). Large plasmid in *Agrobacterium tumefaciens* essential for crown gall-inducing ability. Nature. 252: [169-170](#).

Van Larebeke, N., Genetello, C.H., Schell, J., Schilperoort, R.A., Hermans, A.K., Hernalsteens, J.P. and Van Montagu, M. (1975). Acquisition of tumour-inducing ability by non-oncogenic agrobacteria as a result of plasmid transfer. Nature. 255: [742-743](#).

T-DNA transfer from the *Ti* plasmid into the plant cells (1977)





Figure 7. Preparative Horizontal Agarose Slab Gel Electrophoretic Separation of A277 Plasmid Fragments A277 plasmid was cleaved with *Sma* I


Figure 5. Renaturation Kinetics of *A. tumefaciens* A277 Plasmid *Sma*I Digest Band 3 DNA in the Presence of Tumor DNA and Control DNAs



- Extract tumor DNA
- Hybridize with Ti plasmid



Tumor-inducing and opine synthesis genes on T-DNA imported into plant cells (1977-1985)




T4SS = Type IV Secretion System






Agrobacterium genome components



Agrobacterium tumefaciens C58 strain


Ti plasmid: key components for *Agrobacterium* T-DNA transfer



Jin SG, Prusti RK, Roitsch T, Ankenbauer RG, Nester EW (1990) Phosphorylation of the VirG protein of *Agrobacterium tumefaciens* by the autophosphorylated VirA protein: essential role in biological activity of VirG. *J Bacteriol.* 172(9):4945-50.

Jin SG, Roitsch T, Christie PJ, Nester EW (1990) The regulatory VirG protein specifically binds to a cis-acting regulatory sequence involved in transcriptional activation of *Agrobacterium tumefaciens* virulence genes. *J Bacteriol.* 1990 Feb;172(2):531-7. (identified **vir box**)

Model of the type IV secretion system with the known or suspected localization of its 12 components



Mary C, Fouillen A, Bessette B, Nanci A, Baron C (2018) Interaction via the N terminus of the type IV secretion system (T4SS) protein VirB6 with VirB10 is required for VirB2 and VirB5 incorporation into T-pili and for T4SS function. *J Biol Chem.* 293(35):13415-13426.




Agrobacterium as a vector for plant transformation



The idea of plant transformation was brewing in late 1970's

Schell J, Van Montagu M. (1977) Transfer, maintenance, and expression of bacterial Ti-plasmid DNA in plant cells transformed with *A. tumefaciens*. [Brookhaven Symp Biol.](#) (29):36-49.

Chilton MD. (1979) *Agrobacterium* Ti plasmids as a tool for genetic engineering in plants. [Basic Life Sci.](#) 1979;14:23-31.



Schell J, Van Montagu M, De Beuckeleer M, De Block M, Depicker A, De Wilde M, Engler G, Genetello C, Hernalsteens JP, Holsters M, Seurinck J, Silva B, Van Vliet F, Villarroel R (1979) Interactions and DNA transfer between *Agrobacterium tumefaciens*, the Ti-plasmid and the plant host. [Proc R Soc Lond B Biol Sci.](#) 204(1155):251-66.

Zambryski P, Holsters M, Kruger K, Depicker A, Schell J, Van Montagu M, Goodman HM. (1980) Tumor DNA structure in plant cells transformed by *A. tumefaciens*. [Science.](#) 209(4463):1385-91.

Cloned T-DNA, direct evidence that a specific piece of the Ti plasmid, the T region, can be transferred and integrated into plant cell DNA during tumor formation.

Wullems GJ, Molendijk L, Ooms G, Schilperoort RA. (1981) Retention of tumor markers in F1 progeny plants from in vitro induced octopine and nopaline tumor tissues. [Cell.](#) 24(3):719-27.



Tumorous tobacco shoots
Grafted to non-transformed tobacco




Figure 1 A. Autoradiograph of **Southern Blots** Obtained with EcoRIRI digested Transformant Callus and Shoot DNA and Normal Tobacco DNA Hybridized to a 32P-labeled T-Region Clone

- 1) Ti plasmid
- 2) Control
- 3) Tumor shoot 1
- 4) Tumor shoot 2
- 5) Tumor callus

B

Eco R I	2	6	7	24	13	22	1
Bst-Bam H I	1	8	19	2			

pNW 31C-8,29-1



Figure 3. Mature, Flowering Plants Obtained Following Grafting of Tumorous Shoots and Normal SRI Shoots



Plant transformation by Agrobacterium in 1983

Caplan A, Herrera-Estrella L, Inzé D, Van Haute E, **Van Montagu M**, Schell J, Zambryski P. (1983) [Science](#). 222(4625):815-21.

Partially disarmed pGV3850
pNos-kanR, tobacco protoplast
Co-integration vector

Barton, K.A., Binns, A.N., Matzke, A.J.M. and **Chilton, M.-D.** (1983).
Regeneration of **intact tobacco plants** containing full length copies of
genetically engineered T-DNA, and transmission of T-DNA to R1 progeny.
Cell. 32: [1033-1043](#).

Wide host range plasmid
pRK290:pBR325 **kanR (36 kb)**

Fraley RT, Rogers SG, Horsch RB, Sanders PR, Flick JS, Adams SP,
Bittner ML, Brand LA, Fink CL, Fry JS, Galluppi GR, Goldberg SB,
Hoffmann NL, Woo SC. (1983) Expression of bacterial genes in plant
cells. [Proc Natl Acad Sci U S A](#). 1983 Aug;80(15):4803-7.

Agrobacterium GV3111
Petunia/npt2
co-integration vector



2013 World Food Prize




Monocot transformation by *Agrobacterium* **requires** induction signals

Wang K, Stachel SE, Timmerman B, VAN Montagu M, Zambryski PC (1987) Site-Specific Nick in the T-DNA Border Sequence as a Result of *Agrobacterium* vir Gene Expression. [Science](#). 235(4788):587-91. (**acetosyringone-induced border cleavage**)

Li XQ, **Liu CN**, Ritchie SW, Peng JY, Gelvin SB, Hodges TK. (1992) Factors influencing *Agrobacterium*-mediated transient expression of *gusA* in **rice**. [Plant Mol Biol](#). 20(6):1037-48. (AB salts, 1% glucose, 20 mM MES buffer **pH 5.6**, 2 mM sodium phosphate buffer pH 5.6, 10 mM octopine and **100/μM acetosyringone**)

Chan MT, Chang HH, Ho SL, Tong WF, Yu SM (1993) *Agrobacterium*-mediated production of transgenic **rice** plants expressing a chimeric alpha-amylase promoter/beta-glucuronidase gene. [Plant Mol Biol](#). 22(3):491-506.



(**Co-incubation of potato suspension culture (PSC) with the *Agrobacterium* inoculum**)

Hiei Y, Ohta S, Komari T, Kumashiro T. (1994) Efficient transformation of **rice** (*Oryza sativa* L.) mediated by *Agrobacterium* and sequence analysis of the boundaries of the T-DNA. [Plant J](#). 6(2):271-82.

[Nat Biotechnol](#). 1996 Jun;14(6):745-50. (**maize**)

[Plant Physiol](#). 1997 Nov;115(3):971-980. (**wheat**)

Agrobacterium suspension+**200 μM acetosyringone induction** + immature embryos



Agrobacterium and strain classification

Order: *Rhizobiales*
Family: *Rhizobiaceae*
Genus: *Agrobacterium*

Flores-Félix JD, Menéndez E, Peix A, García-Fraile P, Velázquez E. (2020) History and current taxonomic status of genus *Agrobacterium*. *Syst Appl Microbiol*. 43(1):126046.

Species	Tumor or root inducing	Common strain
<i>Agrobacterium tumefaciens</i>	Ti plasmid	A6, C58, Ach5, T37
<i>Agrobacterium rhizogenes</i>	Ri plasmid	NCPPB2659 (K599), A4
<i>Agrobacterium rubi</i>	Ti	
<i>Agrobacterium vitis</i>	Ti	

Traditionally call *Agrobacterium* strain type based on opine synthesis

A. tumefaciens:

nopaline



octopine

agrocinopine

pTiC58, pTiT37

pTiA6, pTiAch5

pTiBo542



A. rhizogenes (*Rhizobium rhizogenes*, new):



Cucumopine

mikimopine

pRi2659

pRi1724

Disarmed *ABI* strain Ti plasmid



trfA and *kanR* were integrated into GV3101 Ti plasmid (GenBank [NZ_KY000036](https://www.ncbi.nlm.nih.gov/nuclseq/KY000036)) (TrfA is required for maintaining RK2 *oriV* binary vector)



Frequently used disarmed *Agrobacterium* strains


Strain Name	Chromosome background	Ti-plasmid types	Antibiotic Resistance*	Reference
LBA4404	Ach5	a disarmed octopine-type Ti plasmid pAL4404	Rm ^R	Hoekema et al., 1983
EHA101	C58	a disarmed agropine-type Ti plasmid pEHA101 (pTiBo542ΔT-DNA)	Rm ^R , Km ^R	Hood et al., 1986
EHA105	C58	a disarmed agropine-type Ti plasmid pEHA105 (pTiBo542ΔT-DNA)	Rm ^R	Hood et al., 1993
A136	C58	Cured of Ti plasmid	Rm ^R	Watson et al., 1975
AGL0	C58	a disarmed pTiBo542ΔT-DNA	Rm ^R , Cb ^R	Lazo et al., 1991
AGL1	C58	a disarmed pTiBo542ΔT-DNA	Rm ^R , Cb ^R	Lazo et al., 1991
C58C1	C58	Cured of Ti plasmid	Rm ^R	Deblaere et al., 1985
C58C1(pTiB6S3ΔT, pCH32)	C58	a disarmed octopine-type Ti plasmid pTiB6S3ΔT-DNA and a helper plasmid pCH32	Rm ^R , Cb ^R , Tc ^R	McBride and Summerfelt, 1990
GV3101	C58	Cured of Ti plasmid	Rm ^R	Holsters et al., 1980
GV3101(pMP90)**	C58	a disarmed nopaline-type pTiC58ΔT-DNA	Rm ^R , Gm ^R	Koncz and Schell, 1986
GV3101(pMP90RK)**	C58	a disarmed nopaline-type pTiC58ΔT-DNA+trfA	Rm ^R , Gm ^R , Km ^R	Koncz and Schell, 1986
ABI	C58	a disarmed nopaline-type pTiT37ΔT-DNA+trfA	Rm ^R , Gm ^R , Km ^R	Legacy Monsanto

*Cb, carbenicillin; Em, erythromycin; Gm, gentamicin; Km, kanamycin; Rm, rifampicin, Tc: tetracyclin

** Many commercial sources and publications labeled as GV3101 without properly indicating pMP90 helper plasmid. The strain should contain disarmed pTiT37 as we sequenced.

Source: <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC6501860/>



Co-integrate Vector





Binary Vector

Free vector in *Agrobacterium*, replicate in both *E. coli* and *Agrobacterium*



Hoekema, A., Hirsch, P.R., Hooykaas, P.J.J. and Schilperoort, R.A. (1983). A binary plant vector strategy based on separation of *vir*- and T-region of the *Agrobacterium tumefaciens* Ti-plasmid. *Nature*. 303: [179-180](#). [Disarmed *Agrobacterium* LBA4404](#)
Bevan M. (1984) Binary *Agrobacterium* vectors for plant transformation. [Nucleic Acids Res.](#) 12(22):8711-21.

Ye X, Williams EJ, Shen J, Johnson S, Lowe B, Radke S, Strickland S, Esser JA, Petersen MW, Gilbertson LA. (2011) Enhanced production of single copy backbone-free transgenic plants in multiple crop species using binary vectors with a pRi replication origin in *Agrobacterium tumefaciens*. [Transgenic Res.](#) 20(4):773-86.



2T-DNA strategies for marker-free segregation in progeny



Conventional 2T vector



2 right border 2T vector




Tandem 2T vector



2 binary vector

2 *Agrobacterium* co-transformation

Marker gene autoexcision





DISCLAIMER

THE INFORMATION CONTAINED HEREIN IS EXPERIMENTAL IN NATURE AND IS PROVIDED "AS IS". BAYER MAKES NO REPRESENTATIONS AND EXTENDS NO WARRANTIES OF ANY KIND, EITHER EXPRESSED OR IMPLIED, AS TO THE MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, OR THAT THE USE OF THE INFORMATION WILL NOT INFRINGE ANY THIRD-PARTY PATENT, COPYRIGHT, TRADEMARK, OR OTHER PROPRIETARY RIGHTS.



Thank you!



Any questions?



Pioneers in Agrobacterium biology



Marc Van Montagu (left) and Jeff Schell (right), in 1993.



Pat Zmbryski



Rob Schilperoort (1938–2012), Leiden U



Paul Hooykaas



Eugene Nester




Mary-Dell Chilton



Soil bacterium, *Agrobacterium*, a natural genetic engineer

- Reported in 1907 by Erwin Smith and Charles Townsend in *Science*
- Demonstrated bacteria-free tumors at the Rockefeller Institute by Armin Braun in 1940s
- Pathogenicity transferred from a tumor producing pathogenic strain to a non-pathogenic strain, speculated a plasmid but all isolation failed (Allen Kerr Australia 1969)
- Reported 0.2 Mb plasmid by **Jeff Schell and Marc Van Montagu** in Ghent, Belgium & **E. Nester**, M. P. Gordon, and M.-D. Chilton, Seattle in 1974
- T-DNA is flanked by 25 bp direct repeats (Narendra Yadav in Mary-Dell Chilton's group)
- Acetosyringone vir gene induction, ssDNA formation (Scott Stachel, Pat Zambryski)
- Opine genes were transferred in a Mendelian segregation (Léon Otten in Jeff Schell's group)
- Binary vector system with 25 bp border flanks (Mary-Dell Chilton, Rob Schilperoort)
- Transgenic plants in 1983 (the Ghent group, **Rob Fraley**'s group, and Chilton's group)
- Transgenic rice by *Agrobacterium* (Hei et al. Japan Tobacco Co. **1994**)



Recent taxonomy of *Agrobacterium* species



Table 1

Species currently included in the genus *Agrobacterium* and species causing tumours or hairy roots currently included in other genera of Family *Rhizobiaceae*.

Species	Source of isolation	Pathogenicity symptoms	References
Genus <i>Agrobacterium</i>			
<i>A. radiobacter</i> ^a	Soil and plant rhizosphere	Non-pathogenic	[4,9]
<i>A. tumefaciens</i> ^a	<i>Malus</i> sp. tumours	Tumours	[9,46]
<i>A. rubi</i>	<i>Rubus</i> sp. tumours	Tumours	[17,48]
<i>A. larrymoorei</i>	<i>Ficus benjamina</i> tumours	Tumours	[5]
' <i>A. albertimagni</i> '	<i>Potamogeton pectinatus</i>	No data	[42]
' <i>A. fabrum</i> '	<i>Prunus</i> sp., <i>Humulus lupulus</i> , <i>Euonymus alata</i> , <i>Rubus macropetalus</i> tumours	Tumours	[29]
<i>A. pusense</i>	<i>Cicer arietinum</i> rhizosphere	No data	[31,35]
<i>A. nepotum</i>	<i>Prunus</i> , <i>Vitis</i> and <i>Rubus</i> tumours	Tumours	[31,38]
<i>A. skirniwicense</i>	<i>Chrysanthemum</i> and <i>Prunus</i> tumours	Tumours	[39]
<i>A. arsenijevicii</i>	<i>Prunus</i> and <i>Rubus</i> tumours	Tumours	[25]
' <i>A. deltaense</i> '	<i>Sesbania cannabina</i> nodules	No data	[55]
<i>A. salinitolerans</i>	<i>Sesbania cannabina</i> nodules	No data	[56]
' <i>A. bohemicum</i> '	<i>Papaver somniferum</i>	Non-pathogenic	[62]
<i>A. rosae</i>	<i>Rosa x hybrida</i> tumours	Tumours	[26]
Genus <i>Allorhizobium</i>			
<i>A. vitis</i>	<i>Vitis vinifera</i> tumours	Tumours	[31,32]
Genus <i>Rhizobium</i>			
<i>R. rhizogenes</i>	<i>Malus</i> sp.	Hairy roots	[9,40]
' <i>R. tumorigenes</i> '	<i>Rubus</i> sp. tumours	Tumours	[27]

Flores-Félix JD, Menéndez E, Peix A, García-Fraile P, Velázquez E. (2020) History and current taxonomic status of genus *Agrobacterium*. [Syst Appl Microbiol.](https://www.sciencedirect.com/science/article/pii/S0723202019303418?via%3Dihub) 43(1):126046.
<https://www.sciencedirect.com/science/article/pii/S0723202019303418?via%3Dihub>



Protocol training: Agrobacterium

Xudong Ye



The goals of this training

- Prepare *Agrobacterium* glycerol stock for long term storage
- Prepare *Agrobacterium* competent cells
- Transform a binary vector into *Agrobacterium*
- Verify binary vector in *Agrobacterium* (quality control)
- Prepare *Agrobacterium* glycerol stock for plant transformation
- *Agrobacterium* virulence gene induction for transformation



Key points for Agrobacterium handling


- Always culture *Agrobacterium* at 28~32 °C. Ti plasmid may loss at 37 °C
- Induce virulence genes between 19~28 °C. Often 23 °C or room temperature. Some Vir proteins are not stable at 28 °C
- Use minimal salt media for induction. Optimized pH 5.3




Verify and check Agrobacterium strains when you receive a new Agrobacterium culture

Step 1: check antibiotic resistance on solid plates: [GV3101 \(pMP90\)](#): gent^R, spec^S, kan^S
or use specific primers for PCR detection [pTiGV3101](#) (GenBank # [NZ_KY000036](#))



Step 2: Purify strain by single colony



1 loop in 10 ml
medium, gentamicin
30 mg/L, shake for
3~6 hours to get
spread cells



1 ul culture in 10
ml LB dilution, mix
well



Transfer and spread 10 ul
onto LB solid plate with
gent 30 mg/L, culture at
28 °C for 2~3 days for
single colonies



Agrobacterium seed glycerol preparation for storage

Agrobacterium strain **GV3101 (pMP90)**: gent^R, spec^S, kan^S



Pick a single colony
on solid medium

10 ml LB+
gent 30 mg/L

Shake at 28 °C at 225 rpm
for 24-48 hrs, OD₆₀₀+1.0

2.5 ml **80%** glycerol
7.5 ml **Agrobacterium**
Vortex

Aliquot 1 ml in
a cryogenic vial
Store at -80 °C

Prepare *Agrobacterium* competent cells from glycerol stock

Agrobacterium strain **GV3101 (pMP90)**: gent^R, spec^S, kan^S



Transfer a loop
or 50 μ l of
Agrobacterium

50 ml LB+
gent 30 mg/L
Shake at 28C
at 225 rpm
overnight

Transfer 10 or 20 ml
into 200 **LB only**
Shake at 28C
at 225 rpm
until **OD₆₀₀ = ~0.8**

Centrifuge for 5-10
min at ~5000g on ice

Wash pellet with
200 ml Deionized/
sterilized cold
water



Wash pellet with
200 ml cold
10% glycerol



Wash pellet with
200 ml cold
10% glycerol



Suspend in 5 ml
10% glycerol on
Ice. **Aliquot 50 μ l**
store at -80C



Transform a binary vector into *Agrobacterium*

Binary vector antibiotic resistance: *aadA*/spectinomycin


Agrobacterium strain GV3101 (pMP90): gent^R, spec^S, kan^S

Thaw competent cells on ice

Add 1 μ l plasmid (0.1~1 μ g/ μ l) into 50 μ l competent cells




Transfer the mix (100 μ l tip) into 1 mm gap cuvette on ice, covered, insert to the electroporator holder



Set 2 kv for electroporation (1 mm gap cuvette)




Use 1 ml tip, add 1 ml LB to the cuvette and transfer into a labeled tube



Shake for ~2hrs at 28 °C

Transfer 10 μ l and 100 μ l onto 2 plates (LB gent 30, spec 75 mg/L), spread evenly



Incubate at 28 °C for 2~3 days



Three methods to verify binary vector in Agrobacterium cells (QC)


After minipreparation of plasmids from *Agrobacterium* liquid culture

- 1) Submit miniprep plasmid for sequencing (**slow, expensive, most accurate, recommended**)
- 2) Or transfer 1 ul miniprep plasmid into E coli competent cells by electroporation, add 3~5 ml LB with antibiotic selection, shake overnight, miniprep 1.5 ml culture, check with restriction enzymes (**less skill required**)
- 3) **Direct QC Agrobacterium plasmids (more skills)**



Pick 3-8 single colonies into individual tubes

Pick 3-8 single colonies into individual tubes



Sterile a toothpick Or plastic tip




Add 10 ml LB with gent/spec selection shake at 30 °C overnight

Verify binary vector in Agrobacterium

GV3101/pMON-Cpf1 QC by restriction enzymes


Hind3	12.8, 3.6, 1.46, 0.08 kb
NcoI	11.5, 5.3, 1.2 kb
XhoI	7.1, 6.3, 4.2, 0.37 kb



4 ml overnight culture was purified by 300 ul P1/P2/P3 method, dried, suspended in 100 ul water. 10 ul was used for digestion

Mich: Michelle's sample; **1,2,3**: 3 single colonies of GV3101 re-transformation; **P**: ~0.2 ug/lane purified plasmid

Prepare Agrobacterium glycerol stock for plant transformation





Virulence gene induction for maize transformation

1. Grow *Agrobacterium* at 30 °C overnight in 5 ml rich medium (such as LB or YEP) containing the appropriate antibiotics. (gentamicin 30 mg/L for GV3101, spectinomycin 75 mg/L for plasmid)
2. Dilute approximately 0.5 ml of the culture into 50 ml AB-sucrose minimal medium containing the appropriate antibiotics. Grow overnight at 30 °C until the bacteria are in late-log phase (OD600= ~0.8)
3. Spin down the bacteria. Resuspend in two volumes of induction medium containing 100 µM acetosyringone. Shake 14-24 hours (overnight) **at room temperature** (not at 30°C).
4. Spin down the bacteria. Resuspend in MS plant tissue culture medium. Inoculate plants.
5. After two days, rinse the plant tissue in medium containing 100 µg/ml timentin. Continue incubating the tissue on solidified medium containing timentin.
6. Stain the tissues after various periods of time (2-10 days) in X-gluc.



DISCLAIMER

THE INFORMATION CONTAINED HEREIN IS EXPERIMENTAL IN NATURE AND IS PROVIDED "AS IS". BAYER MAKES NO REPRESENTATIONS AND EXTENDS NO WARRANTIES OF ANY KIND, EITHER EXPRESSED OR IMPLIED, AS TO THE MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, OR THAT THE USE OF THE INFORMATION WILL NOT INFRINGE ANY THIRD-PARTY PATENT, COPYRIGHT, TRADEMARK, OR OTHER PROPRIETARY RIGHTS.



Bayer Russia Plant Biotechnology Conference:

Day 1	Plant Biotechnology Overview and Evolution
Day 2	Introduction to Plant Vectors & Agrobacterium
Day 3	Introduction to Plant Transformation - Corn
Day 4	Soy Plant Transformation & Gene Expression
Day 5	Introduction to Plant Health & Controlled Environment (CE)



Thank you!



Any questions?

