



Plant Biotechnology

Bayer Russia Plant Biotechnology
Conference

July 2023





Molecular Assays

Edit Detection



Gene Editing Detection and Analysis

Starting Materials:

1. Sequencing Data
2. Target sequences
3. Query sequences
4. Metadata



Desired Output:

1. A report of editing events
2. A report on quality of input data
3. Provenance



Corn Zm 7.1 target information

gRNA target highlighted

Gene-specific primer underlined

>Zm7.1_amplicon

CGAAGTACCATGTTTCAAATGCTGGCCGGCTACTGCCTAA**GAGGGCATGCCATCATATTATAC**TAAATAGACTATCATAAAGAACAGATGGTATTATTTC
ACATCGAATGACTCATTGTCCTCGCTAGCACTGTCAACATAAACCAATGTTATTAGAGCCAATTTCCGCGCTGCTGAAGAATGAAGACTGAGCTTGCA
CATCCAGAAGCATCCATCACATGCTTTGTTTTCTTCTTCTTACATTTTCAAGTCCTCCCAACTCACTCCC



Soy Fad2-1A target information

gRNA target highlighted

Gene-specific primer underlined

>Expected Chromosome 10 Fad2-1A_amplicon sequence from gene-specific primers

CCATGCCTTCAGCAAGTACCA**A**ATGGGTTGATGATGTT**G**TGGGTTTGACCC**CTT**CACTCA**A**CACTTTTA**GTCCCTTATTTCTCATGGAAA**
AT**A**AAGCCATCGCCGCCA**T**CACTCCAACAC**A**GGTTCCCTTGACCGTGATGAAGTGTTTGTCCCAAACCAA**A**ATCCAAAGTTGCATG
GT**TTT**CCAAGTAC**TTA**AACAACCCTCTAGGAAGGGCTG**TTT**CTCTTCTC**G**TCACACTCACAATAGGGTGGCCT**AT**GTATTTAGCCTT
CAATGTCT


>Expected Chromosome 20 Fad2-1A_amplicon sequence from gene-specific primers

CCATGCCTTCAGCAAGTACCA**C**ATGGGTTGATGATGTT**A**TGGGTTTGACCC**G**TTCACTCA**G**CACTTTTA**GTCCCTTATTTCTCATGGAAA**
AT**A**AAGCCATCGCCGCCA**C**CACTCCAACAC**G**GGTTCCCTTGACCGTGATGAAGTGTTTGTCCCAAACCAA**A**ATCCAAAGTTGCATG
GT**ACA**CCAAGTAC**CTG**AACAACCCTCTAGGAAGGGCTG**CTT**CTCTTCTC**A**TCACACTCACAATAGGGTGGCCT**TT**GTATTTAGCCTT
CAATGTCT

- There are 14 SNPs that can be used to differentiate edits between Chr10 and Chr20 (SNPs bold and italicized)



1. Perform quality control on sequencing data



- We want to exclude the portions of the reads where quality drops
- Quality of editing calls is directly related to quality of reads being used


F = forward
R = reverse



2. Align sequencing data to target

Reference

AGTCGCATAGCGAGTCGCATAGCGCATCGAGCATAGCGCATCGATTAGGGATATCTAGAGACTTTCTGAAGC



F = forward
R = reverse



3. Find and count edits

For each read see if there is an edit:

AGTCGCATAGCGAGTCGCATAGCGCATCGAGCATAGCGCATCGATTTAGGGATATC

Read123-F

AATCAAGTCGCATAGCGAGTCG---

CGAGTCG-----CGCATCGAGC

Read123-R

+1 count for this deletion

F = forward
R = reverse



4. Report results

Sample	Target	Query	Type	Edit	%	Sequence
Sample001	Zeamays_chr3	Target1	Wt	-	99%	ACGTACGTA ACGTACGTA
Sample001	Zeamays_chr3	Target1	Edit	3M5D1M	15%	ACGTACGTA ACG-----A
Sample002	Zeamays_chr4	Target2	Edit	4M1D4M	40%	GTAGTCTTT GTAG-CTTT

A string annotation for the edit where
M is “match” and D is “deletion”

Alignment to reference



5. Continue with promising edits

- At this point you can decide which edits are promising and move them forward in the pipeline
- Determine which editing technology is working well or not for specific targets
- Add this data to previous experiments for continued provenance tracking of an edit through generations



Available Applications



CRISPResso2

Analysis of genome editing outcomes
from deep sequencing data

- <http://crispresso2.pinellolab.org/submission>

Cas-Analyzer


- <http://www.rgenome.net/cas-analyzer>

CRISPRmatch

- <https://github.com/zhangtaolab/CRISPRMatch>

Hi-TOM 2.0

- http://hi-tom.net/#/mutation/mutation_detection





DISCLAIMER

THE INFORMATION CONTAINED HEREIN IS EXPERIMENTAL IN NATURE AND IS PROVIDED "AS IS". BAYER MAKES NO REPRESENTATIONS AND EXTENDS NO WARRANTIES OF ANY KIND, EITHER EXPRESSED OR IMPLIED, AS TO THE MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, OR THAT THE USE OF THE INFORMATION WILL NOT INFRINGE ANY THIRD-PARTY PATENT, COPYRIGHT, TRADEMARK, OR OTHER PROPRIETARY RIGHTS.



Thank you!



Any questions?

